## REF XS50

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## ENGLISH INSTRUCTIONS FOR USE

### A. DESCRIPTION

The SINAPI Chest Drain is a chest drainage unit incorporating a dry seal, small fluid collection reservoir, air leak detector and suction bulb.

### B. INDICATIONS

The SINAPI Chest Drain is used for:

- The evacuation of air from the pleural cavity.
- The prevention of air re-accumulation in the pleural cavity.
- The facilitation of lung re-expansion and restoration of normal breathing dynamics.

### C. SETUP

HOW TO CONFIRM CHEST CATHETER POSITION DIRECTLY AFTER INSERTION

- Ensure a light fit between the inlet tubing connector and the chest catheter.

To confirm **Pneumothorax**: Look for vapour on the inside of the chest tube. Depress the bulb, the bulb will not stay depressed but initially re-expand. If it does not re-expand it may confirm an absence of air at the catheter tip. This will mean that the catheter is **not placed in the pleural space** (or pneumothorax is resolved).

**MOBILE:** The short connecting tube facilitates mobilisation of patients.

### D. MANUALLY INCREASING THE DRAINAGE RATE

Increase the drainage rate by depressing the bulb until the bulb stays depressed Fig. (1). If the bulb re-expands, pinch the tubing above the valve and then depress the bulb. Release both. Repeat. This applies suction and increases the drainage rate.

For BP fistula (persistent air-leak) pneumothorax will remain and the bulb will not stay depressed - attach external suction if required by physician.

### E. ATTACHING TO WALL SUCTION

If the bulb does not stay depressed, attach suction. Set suction between -5cm H2O and -45cm H2O pressure - as prescribed by physician.

### F. DETECTING BUBBLING / PERSISTENT AIR-LEAK

Depress the bulb. If it stays depressed there is no air leak. If it re-expands, syringe 20ml into the air leak chamber via the needle-free port.

**Bubbling or saline level movement = air leak** Fig. (2) (Air-Leak Chamber)

To monitor air leak, ensure that the saline in the Air Leak Chamber is maintained to the fill line at all times.

### G. MONITORING

Fig. (3a) (Good News) DEPRESSED BULB - Reason: Suction/Negative pressure is achieved - provided there is no blockage in tubing

Fig. (3b) (Bad News) EXPANDED BULB - Potential Reasons: 1) Pneumothorax 2) Fistula (Persistent Air Leak) 3) Tube connections are not secure (air leaks into the system)

### H. A GUIDE WHEN TO REMOVE THE CATHETER (You need to confirm both 1. and 2.)

- Drainage volume must be low.
- Depress the bulb firmly: the bulb must stay depressed and the valve angled (indicating negative intra-pleural / mediastinal pressure). **Confirm clinically and/or with chest x-ray.**

### I. MAINTENANCE

- A small volume of liquid might accumulate in the tube above the valve. This is normal and an indication of negative pressure above the valve.
- During outpatient management, secure the chest drain reservoir to the patient to prevent dislodgement of the catheter.

### K. CAUTIONS and WARNINGS

- This device is only intended for use by appropriately trained personnel.**
- Compatibility of the Sinap Chest Drain with thoracic catheters needs to be established by the user.
- Assess the thoracic catheter connections for signs of airn air leak.
- Thoracic catheter dressings should be occlusive.
- Ensure that there are no dependent loops or kinks in the patient drainage tube.
- The volume, rate and nature of the fluid collected in the SINAPI Chest Drain should be monitored regularly.
- Little or no drainage may indicate a blockage in the drainage system.
- Do not clamp tube, this will inhibit drain operation and may compromise respiratory function of patient.
- Samles of the drained fluid can be taken directly from the needle-free sampling port located on the front of the device.
- Ensure that there are no dependent loops or kinks in the patient drainage tube.
- The volume, rate and nature of the fluid collected in the SINAPI Chest Drain should be monitored regularly.
- Connect only ONE chest catheter per chest drain.

## PATIENT INFORMATION FOR HOME CARE

### A. EMPTYING THE DEVICE

- Wash your hands for at least 30 seconds & dry with a clean towel.
- Wipe the Needle-free sampling port with an alcohol swab & let it dry.
- Gently push and twist the tip of a new syringe into the port.

- Titl the device to the left to allow fluid to pool at the sampling port.
- Pull the plunger back on the syringe to remove the fluid.
- Remove the syringe & empty the drainage into your toilet.
- Discard the empty syringe into the bin.

- Wash and dry your hands.
- Record the colour and amount of the fluid that you removed.

### IMPORTANT

- If the fluid level in the collection chamber exceeds the Max fill line, it may leak out.
- Empty the device when it reaches the Max fill line and preferably before you lie down
- To prevent infection, do not touch the syringe tip or the port when emptying the device
- Your physician or nurse will inform you regarding dressing changes.

### B. SHOWERING AND BATHING

The device should not be submerged - it is best to either wash at a sink or take a shower. When showering, cover the dressing with plastic and secure it with tape to keep it dry.

### C. INFORM YOUR PHYSICIAN & GO TO THE CLOSEST EMERGENCY UNIT IF:

- The chest tube accidentally comes out:
  - Do not try to reinsert the tube.
  - Immediately cover the chest wound with a dressing and plaster.
- The device disconnects from the chest catheter:
  - Knk it and try to reconnect.
  - Press the bulb a few times after reconnecting it.
  - Keep the tube kniked if you are unable to reconnect it.
- You experience any **breathing difficulty** or a **new chest pain**.
- Fluid drainage **increases / stops** or **changes colour**.

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